IVD

CE

nova Max Pro™ CREAT
eGFR

**Creatinine and eGFR Meter System
Quick Reference Guide**

nova[®]
biomedical

Quick Reference Guide	1
Kurzanleitung	17
Οδηγός γρήγορης αναφοράς	33
Guia de referencia rapida	49
Guide de référence rapide	65
Guida rapida	81
Korte handleiding	97
Guia rápido de referências	113
Gyors referencia kézikönyv	129

Refer to the Nova Max Pro™ Creatinine and eGFR Meter System Instructions for Use Manual and Package Inserts for complete instructions for use, indications, precautions, and limitations of the system.

Intended Use

The Nova Max Pro Creatinine and eGFR Meter System is intended for in vitro diagnostic use for the quantitative measurement of creatinine in capillary whole blood obtained from the fingertip for screening of kidney health in asymptomatic persons with potential kidney disorders such as Chronic Kidney Disease (CKD) and Acute Kidney Injury (AKI) by calculating Estimated Glomerular Filtration Rate (eGFR). It is intended for multiple patient use in Near-Patient/Point of Care and professional healthcare settings. The system should only be used with the Nova Max Pro Creat eGFR Creatinine Test Strips and Nova Max Pro Creat eGFR Creatinine Control Solutions and with single-use, auto-disabling lancing devices. It is not intended to diagnose a specific kidney disease or condition. The system should not be used to alter kidney disease treatment by changing any medication schedule or dosage unless specifically instructed by a healthcare professional.

WARNING: Blood samples and blood products are potential sources of hepatitis and other infectious agents. Handle all blood products with care. Wear gloves when performing measurements on another person. Items that are used to measure Creat and eGFR, i.e., Test Strips, lancets, and alcohol swabs, must be disposed of in accordance with local regulations to avoid risk to anyone.

Limitations

The Nova Max Pro Creat eGFR Creatinine Test Strips give accurate results when the following limitations are observed:

- Blood – Use only capillary whole blood. Do not use serum or plasma.
- Do not use the Nova Max Pro Creatinine and eGFR Meter System for testing neonates.
- Test results are best obtained when used within an operating relative humidity of 10% to 90% (non-condensing) and an operating temperature range from 15 to 40°C. Testing outside this range may cause inaccurate results.
- Altitude – There is no effect of altitudes up to 12,000 feet (3658 meters) above sea level.

General Safety

Personnel operating this Meter must be proficient in the operating and maintenance procedures of the Meter, the following safety procedures must be followed.

1. Read the safety and operating instructions before operating the Meter.
2. Retain the safety and operating instructions for future reference.
3. Observe all warnings on the Meter and in the operating instructions.
4. Follow all operating and use instructions.

Testing a Blood Sample


1. From the Home screen, press Run Blood Test.
2. The Enter Patient ID screen displays. You can enter the ID by number only.

3. Enter Patient ID allows the patient's ID to be stored with their data.
4. Press the Continue button to display the Enter Patient Age screen. Enter the patient's age for this test from 18-120 years.
5. Press the Continue button to display the Select Gender screen: press Male or Female.

6. Press the Continue button to display the Select Ancestry screen: press Black or Other.

NOTE: The select Ancestry screen is only displayed if the meter is set to report eGFR results using the CKD EPI (2009) equation.

7. Press the Continue button to display the Insert Strip screen. The meter waits for you to insert a strip. Once the test strip is inserted the meter transitions to the Apply Sample screen.

8. Wash patient's hand with water then dry thoroughly. Alternatively, use alcohol pads to clean area; dry thoroughly after cleaning.
9. Holding hand downward, massage finger with thumb toward the tip to stimulate blood flow.
10. Use a single-use, disposable safety Lancet to puncture the finger. To reduce the risk of pre-analytical error, facilities should consider using a 21-gauge lancet when collecting capillary samples..
11. Squeeze the finger to form a drop of blood. Wipe away the first drop of blood then squeeze the finger again to form a second drop of blood.

12. Touch the end of the Test Strip to the fresh blood drop until the Test Strip fills and the on-screen countdown timer begins. (Beeper sounds if enabled.)
13. A countdown on the screen appears while the test is in progress. The eGFR and Creat results are available on-screen in 30 seconds.
14. The Result Screen displays the eGFR and Creat results. At this point, the test is finished. To save the result to memory, press the Accept button. To delete the result, press the Reject button. The result can be viewed later in the Result screens. You can also click on the personal settings icon  which shows their settings for eGFR. Once you press Reject or Accept, the meter displays the Home screen.
15. The result is automatically stored in memory.

16. Wash your hands thoroughly with soap and water after handling the Meter, Lancet, or Test Strips.

NOTE: Do not press the Test Strip directly against the skin. Touch the Test Strip gently to the blood drop.

NOTE: Discard used Lancets and Test Strips in a puncture-proof container such as a biohazard container.

Red <30 mL/min/1.73m²	Orange 30-44 mL/min/1.73m²	Yellow 45-59 mL/min/1.73m²	Green 60-90 mL/min/1.73m²
---	---	---	--

NOTE: The eGFR number indicates how well the kidneys are functioning. As kidney disease progresses, eGFR goes down.

Testing a Control Solution

NOTE: Read the package insert sheets for Nova Max Pro Creat eGFR Control Solution for complete instructions, indications, precautions, and limitations of the system.

1. From the Home screen, press Run QC Test. The Select QC screen displays. Select Run QC 1, Run QC 2, or QC 3.
2. Once the QC test is selected, the Insert Test Strip screen displays. Insert a test strip.
3. The Add QC 1, QC 2, or QC 3 screen displays. Add the QC solution as shown on the display to the strip.

4. A countdown timer showing the time remaining until test results will be ready is displayed. Once the countdown reaches 0, the meter displays the QC Results screen.
5. The QC Result Screen displays the QC test result. Press the Home button to return to the Home screen. QC results can be recalled if desired from the QC History screen.

NOTE: *If the strip is removed before the test starts or is not used for over 2 minutes, the screen will go blank.*

NOTE: *A quick beep sounds when sufficient Control Solution has been added to the Test Strip.*

NOTE: *Nova Max Pro Creat eGFR Control Solution test results are available on-screen in 30 seconds.*

Out-of-range results may be caused by the following:

- You may not be doing the test properly. Retest and follow the instructions carefully.
- The control solution may have expired or have been contaminated. Check the expiration date on the control solution vial. The control solution is good for only 3 months after opening. Make sure the control solution vial is closed when not in use.
- The test strip may have expired. Check the expiration date on the test strip vial.
- The test strip may have been damaged. Retest using a new test strip.
- The Nova Max Pro Creat eGFR Meter may not be working properly.
- Acceptable quality control ranges may differ at your institution from the ranges reported by Nova Biomedical. Follow your institution's guidelines for quality control results.

NOTE: *If the Control Solution test result is outside the range (is either higher or lower), the Meter and Test Strip may not be working as a system. Repeat the process using a new Test Strip. Do not use the Meter until test results fall within the appropriate range. If the problem continues, call Customer Service.*

Cleaning and Disinfecting the Meter

WARNING: *Cleaning is not the same as disinfecting. Cleaning means removing protein or other contaminants from the surface. Disinfecting means to kill or prevent the growth of disease-carrying microorganisms.*

The Nova Pro Creat eGFR Meter should be cleaned and disinfected weekly over the intended use-life of the Meter. The cleaning and the disinfecting procedure was validated a total of 10,950 times by Nova Biomedical. The validation testing corresponded to cleaning and disinfecting 10 times a day for 3 years.

Meter Cleaning and Disinfection Procedure

The Nova Pro Creat eGFR Meter should be cleaned and disinfected for each patient use to minimize the risk of transmission of blood-borne pathogens between the patient and HealthCare Professionals.

HealthCare Professionals and others should follow Good Laboratory Practice guidelines and these important safety instructions. HealthCare Professionals should ensure that they are wearing protective gloves when disinfecting the meter and should wash their hands thoroughly with soap and water after handling the meter

NOTE: Nova recommends cleaning and disinfecting the Meter with the following EPA Registered products - Clorox® Germicidal Wipes, EPA Registration #67619-12. Super Sani-Cloth® Germicidal Disposable Wipes, EPA Registration #9480-4. These products are widely available at numerous retail and wholesale providers .

NOTE: Cleaning and disinfection may in rare cases damage the device(s). Meter damage may include plastic housing cracks, cloudiness, or frosting of the display, legibility or response issues with the keypad, or battery compartment fluid leakage. Signs of Meter performance deterioration may include failure to recover proper control results or the inability to perform a blood test. If you observe damage due to cleaning and disinfecting, please stop using the Meter and contact Customer Service.

Cleaning and Disinfecting the Meter

NOTE: *To properly clean and disinfect the Meter, steps 1 to 5 should be performed together before testing on each patient.*

1. Clean the Meter.

- Wipe the external surface of the Meter thoroughly with a fresh Germicidal Wipe.
- Discard the used wipe per Step 4.

2. Disinfect the Meter.

- Remove another fresh Germicidal Wipe from the canister. Thoroughly wipe the top, bottom, left, and right sides of the Meter avoiding the Test Strip port by wiping the surface a minimum of 3 times horizontally followed by 3 times vertically.

3. Observe surface contact time.

- Ensure the Meter surface stays wet **for the recommended time** and is allowed to air dry for an additional **1 minute** (1 minute for Clorox Germicidal Wipes or 2 minutes for Super Sani-Cloths).

NOTE: *If you must rewet the surface of the Meter, use a new, fresh wipe.*

4. Dispose of wipes

- Dispose of the used germicidal wipes in a standard waste container.

5. Wash and sanitize hands.

- Wash your hands thoroughly with soap

Error Codes

E2 Operating Temperature Error: Error means no blood test result. Move the meter to an area where the temperature is acceptable (59°F to 104°F or 15°C to 40°C), allow meter to adjust to the temperature. Repeat the test.

E3 Used Strip Error: Test strip is defective or used. Repeat the test with a new test strip. If the error persists, perform the test using an alternate test strip vial or alternate method.

E4 Short Sample Error: An insufficient sample (control or blood) has been placed on the test strip. Repeat the test with a new test strip. If the error persists, perform the test using an alternate test strip vial or alternate method.

Technical Assistance

For technical assistance inside of the U.S., call Nova Biomedical Technical Services at

Phone.: 1-800-545-NOVA 1-781-894-0800

FAX: 1-781-894-0585 Web: www.novabiomedical.com

For technical assistance outside of the U.S., call your local Nova subsidiary or authorized distributor.

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Nova Biomedical GmbH Hessenring 13A, Geb. G
Nova Biomedical Corporation
200 Prospect Street Waltham, MA 02454 U.S.A.